



ATMOS MedizinTechnik GmbH & Co. KG

Ludwig-Kegel-Str. 16
79853 Lenzkirch / Deutschland

Phone +49 (0)7653-689-222
E-mail service@atmosmed.de

PR2GB.110150.0 2022-11 Index: 01 Subject to alterations; errors excepted. Corian® is a registered trademark of E. I. du Pont de Nemours and Company or its affiliated companies.



CORIAN®

The cleaning and care of CORIAN®

*Exactly
what
I need*

OVERVIEW — CARE INSTRUCTIONS

One surface. Many benefits.

Please observe the following instructions for cleaning your CORIAN® worktop or other surfaces made of CORIAN®.

DAILY CLEANING

For everyday cleaning, all liquid residues should be removed with a soft cloth. Alternatively, you can use a damp cloth with a mild abrasive liquid cleanser. Then rinse the surface with warm water. Finally, use a soft cloth to wipe the surface dry.

CORIAN® is a homogeneous, durable, and repairable material. Due to these properties, it can be easily renewed. Generally, most damage such as that caused by severe impact, heat, or chemicals can be repaired on-site. As a result, the surface regains its original, smooth and hygienic quality.

⚠ WARNING

- Please make sure to wipe up spilled liquids immediately. The surface should be cleaned in a circular motion.
- Always use a hot pad or trivet (with rubber feet) to protect the surface from hot objects. Do not place hot objects directly on the CORIAN® surface. Heat exposure can damage the surface.
- Do not pour boiling liquids into the sink or lavatory without running cold water at the same time.
- Please note that under normal daily use, CORIAN® may show slight signs of wear after a certain period of time.
- Do not use any pointed or sharp objects directly on the worktop. Always use a cutting board.
- Dark, heavily pigmented colors are more likely to show scratches, dust, and normal wear and tear than light, textured colors (please refer to our CORIAN® Solid Surface portfolio).
- Make sure that accidental chemical spills are rinsed off immediately using plenty of soapy water. Unnoticed or prolonged contact with chemicals can damage the surface.

ABCs of stain removal.

The following steps apply for all stains:

1. Remove residue with a soft cloth.
2. Use a damp cloth and a mild abrasive liquid cleanser.
3. Rinse the surface several times with warm water and wipe dry with a soft cloth.

The list below contains additional cleaning methods for certain stains.

VINEGAR, COFFEE, TEA, LEMON JUICE, DYE

Step 2. Non-scratch scrub sponge and ammonia-based cleaner or solid surface cleaner.
Alternatively: non-scratch scrub sponge and diluted household bleach.

GREASE, OIL RESIDUES

Step 2. Non-scratch scrub sponge and ammonia-based cleaner or solid surface cleaner.

LIME, SOAP, MINERAL DEPOSITS

Step 2. Rub over the stain with a non-scratch scrub sponge and standard household limescale remover or vinegar.

LIGHT SCRATCHES, MINOR BURN MARKS

Step 2. Non-scratch scrub sponge and diluted household bleach.

INK, MARKING PEN, BLOOD

Step 2. Non-scratch scrub sponge and ammonia-based cleaner or solid surface cleaner.
Alternatively: non-scratch scrub sponge and diluted household bleach.

IRON, RUST

Step 2. Non-scratch scrub sponge and metal cleaner or rust remover.

IODINE, MOLD, MILDEW

Step 2. Non-scratch scrub sponge and diluted household bleach.

If stains cannot be removed or if you are interested in having your CORIAN® work surface sanded down professionally, please contact us by telephone at 07653-689-222 or by e-mail at service@atmosmed.de.