


Quick Guide ATMOS® C 341 Battery

i This Quick Guide does not replace the operating instructions.


For your safety

- Recharge the battery at the latest after 3 months, even if you do not use the device.
- Perform a function check after each use. Perform a function check every 4 weeks in case you do not use the device for a longer period.
- The device may only be used by persons who were medically trained, and were trained in the medical suction.
- Do not operate the device without a bacterial filter.

Device overview




- 1** Control panel
- 2** Battery compartment cover
- 3** Device base
- 4** Release button wall and device support
- 5** Hose rewind with suction hose
- 6** Connection charging accessories
- 7** Guide for the wall and device support



- 8** Display of battery status
- 9** Vacuum controller
- 10** Vacuum display
- 11** On/off button
- 12** Power indicator



Display of the battery status

- The following display values are not valid during battery charging.
- The battery status must be supervised visually.
- If only the red LED lights up, although the battery has been loaded for a longer period, the battery is defective.




Support

Mounting



Removal



Switching on and operate

Switching on



1. Press the on/off button.
 - » The pump starts.
 - » The LED above the on/off button is illuminated as long as the device is switched on.

Vacuum adjustment



2. Turn the vacuum controller for adjusting the required vacuum.
 - » The vacuum display shows the current suction performance.

Switching off



3. Press the on/off button for at least 1 second.
4. Clean and disinfect the device after every use.
5. Perform a function check after each cleaning and disinfection.

Cleaning and disinfection

What	How			When		
	C	D	S	Daily*	Every 14 days*	After each patient / after each suction
Parts to be reprocessed	Cleaning	Disinfection	Sterilization	Daily*	Every 14 days*	After each patient / after each suction
Surfaces						
Housing	X	X		X		X
Device base	X	X ¹		X		X
Hose rewind	X	X ¹		X		X
Wall and device support	X	X				X
Power supply and recharging unit	X	X				X
Canister system						
DDS secretion canister	X	X ¹	X ²	X		X
Outer canister lid (DDS)	X	X ¹	X ²	X		X
Inner canister lid (DDS)	X	X ¹	X ²	X		X
Float ball (DDS)	X	X ¹	X ²	X		X
Filter holder	X	X ¹	X ²	X		X
Sealing ring	X	X ¹		X		X
Bacterial filter					X	X
Fingertip				X		X
Support for canister system	X	X ¹		X		X
Suction hose (DDS)	X	X ¹	X ²	X		X
	X	X ¹	X ²	X		X

1) Preferred: machine cleaning and disinfection in the washer disinfector (rinse cycle up to 95°C)

2) If required hot-steam sterilization at 134° C with devices acc. to EN285

Function check

Manual

1. Please check whether the following parts are damaged or torn:
 - all hoses
 - Canister system
2. Switch on the device.
3. Check: LEDs and battery status.
4. Connect a fingertip to the suction hose and close the auxiliary air vent.
5. Close the front opening of the fingertip with your thumb.
6. Select the vacuum -0.5 bar.
7. Please check whether the device reaches the vacuum after approx. 20 seconds: The pump switches off and the vacuum display shows -0.5 bar.
8. In case the device does not reach the vacuum within 20 seconds: Check the device on potential sources of error.