Authorized distributor / Customer services

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	<b>MANUFACTURER</b> Ulrich GmbH & Co. KG Bachstr. 109/110 89081 Ulm - Germany	<b>DISTRIBUTED IN THE USA BY:</b> Ulrich Medical USA, Inc. Chesapeake, MD 20765		
<b>REF UT 1330-XL CONICAL</b>				

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## EN Disposable cuffs

### Intended use

The intended use of the single use tourniquet cuffs is to regulate the flow of arterial blood in the upper and lower extremities to the point of arresting the flow entirely. This is accomplished by connecting the single use tourniquet cuff to a compatible tourniquet, which inflates the cuff with air. The level and duration of pressure application is up to the sole discretion of the physician. The cuff is used on extremities in neurosurgery as well as for applications in vascular, trauma and plastic surgery.

### Cuff use

- Before use: Make sure that the disposable cuffs are in good working order.

Before arresting the blood flow, the physician must first assess the patient's medical history. Proper selection and use of the tourniquet cuff as well as the level and duration of the pressure applied all depend on the patient's condition (condition of the skin, tissues and vascular system, blood pressure, size and shape of the extremity, general condition, obesity, etc.) and are up to the sole discretion of the physician. The cuff may not work reliably if the end of the cuff does not reach at least to the beginning of the loop side of the Velcro strip.

The disposable cuff is delivered in sterile packaging (ETO) and may only be used once – the cuff must not be reused. Do not use the cuff if the packaging has been damaged or if the sterility date has expired.

### ⚠ Warnings

- The cuff must only be applied by professionals and/ or under the supervision of professionals.
- The various cuff sizes are color coded. The UT 1315 measuring tape is recommended as a tool for determining the proper size.

- Insert thin padding (cotton or other similar material) under the cuff to prevent chemical burns and skin irritation caused by preoperative skin treatment agents or rinse solutions used during surgery. The cuff may also need to be sealed off distally.
- Make sure that no folds of skin form beneath the cuff. This means that the cuff should not be twisted or slid once it is in place and inflated.
- The cuff should be placed in the most proximal position possible, as this will maximize the amount of muscle and tissue mass available between nerves and vascular structures and the cuff.
- In order to monitor the duration of the blood arrest, the cuff should only be used in conjunction with a tourniquet that has been equipped both with safety and alarm features, as well as with a built-in timer.
- Do not use the cuff if the packaging has been damaged or if the sterility date has expired.
- The cuff is designed to be used only once and must not be processed for multiple use. Processing and sterilizing the cuff would destroy its ability to function properly.

### ⚠ Warnings (specific to the XL conical disposable cuff)

- The XL conical disposable cuff has been designed for conical extremities in obese and muscular patients, as well as in patients with water retention in their extremities (lipedema).
- The XL conical disposable cuff must only be placed on the legs.

<p><b>MANUFACTURER</b> Ulrich GmbH &amp; Co. KG Suchbaurweg 12 82081 Ulm   Germany</p> <p><b>DISTRIBUTED BY THE USA BY:</b> Ulrich medical USA, Inc. Chesapeake, MD 20705</p> <p><b>REF</b> UT 1330-XL <b>CONICAL</b></p> <p>CE 0123    </p>	<p><b>Positioning the conical disposable cuff</b></p> <p>This additional symbol helps users position the cuff by depicting how the end of the cuff with the smaller circumference must be placed at the distal end of the extremity. This is the only orientation in which the conical shape of the cuff can be ideally fit to the shape of conical extremities.</p>
--	--

Do not re-use	Latex free	Consult instructions for use	Do not use if package is damaged	Sterilized using ethylene oxide	-	Manufacturer

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