

ATMOS LS Battery LED

Battery-operated light source for endoscopy

*Exactly
what
I need*



Light, small & powerful.

Maximum mobility with the ATMOS LS Battery LED.

- High-performance LED (3 levels of brightness + power mode)
- Long battery life
- Light and easy to handle
- Waterproof
- Magnetic USB charging chip with display of battery status



TECHNICAL DATA

ATMOS LS BATTERY LED 507.4670.0

| | |
|---------------|--------------------|
| Size | 75 mm x 21 mm |
| Weight | 51 g incl. battery |
| Power | 10 watts |
| Charging time | 3 hours |

LEVEL OF BRIGHTNESS

| Power mode | Level 1 | Level 2 | Level 3 |
|-------------|------------|------------|-----------|
| 160,000 lux | 88,000 lux | 12,000 lux | 1,700 lux |



MOBILE ENDOSCOPY

The ATMOS iQam (USB camera) and ATMOS LS Battery LED complement each other perfectly to form a mobile endoscopy set.

