

Flexible
and
safe


High
level of
hygienic
comfort

Equipment
rails for every
requirement

ATMOS Equipment rail

**Secure
mounting for
your devices**


*Exactly
what
I need*



ATMOS Equipment rail

FLEXIBLE AND SAFE

The ATMOS Equipment rails provide you with a flexible solution for securely mounting your medical devices and accessories. Our equipment rails are standardised according to DIN EN ISO 19054 and have a profile of 25 x 10 mm.



HIGH LEVEL OF HYGIENIC COMFORT

The enclosed cover caps not only ensure that the potential equalisation is concealed optically, but also provide a high level of hygienic comfort in combination with the preinstalled end caps.

EQUIPMENT RAILS FOR EVERY REQUIREMENT

Individuality is very important to us. In addition to the standard lengths, we also offer you the option of configuring the ATMOS Equipment rails according to your requirements and wishes. Not just the lengths, but also the distances between the individual wall holders are variable and can be optimally adapted to your needs.

Fixed length

Fixed length	Article number
Equipment rail 25x10 mm / 1 m	401.0702.0
Equipment rail 25x10 mm / 2 m	401.0703.0
Equipment rail 25x10 mm / 3 m	401.0704.0

Variable length

Variable length	Article number
Equipment rail 25x10 mm / 0–1.00 m	401.0705.0
Equipment rail 25x10 mm / 1.01–1.50 m	401.0706.0
Equipment rail 25x10 mm / 1.51–2.00 m	401.0707.0



MedizinTechnik

Our ATMOS
product specialists
are available for
individual
consultation.

ATMOS MedizinTechnik GmbH & Co. KG

Ludwig-Kegel-Str. 16
79853 Lenzkirch / Germany

Phone +49 (0)7653-689-0
Fax +49 (0)7653-689-190
E-mail info@atmosmed.com

www.atmosmed.com

Your local contact: