

ATMOS Hygrotherm plus

# When precision counts!



## ATMOS Hygrotherm plus

ATMOS Hygrotherm plus ear irrigation		REF
<p><b>1 Irrigation unit for the ENT field</b>  with electronically controlled flow heater, irrigation mode (37°C, up to 450 ml/min), energy-saving mode to reduce energy consumption, highest temperature accuracy and safety for the patient. Water supply hose G3/4", 3 m, waste water hose G1/4", 3 m, 220-240 V~, 50/60 Hz.  <i>Included in delivery:</i>  3 x nozzle attachment -II straight, 80 mm; 1 x Plexiglas splash protection;  1 x seal for splash protection; 1 package of hose tips for nozzle (30 pcs.);  1 x water filter housing R.V.-A1/4"; 1 x water filter CellPor-Pe; 1 x mains cable.</p>		502.0901.0
ATMOS Hygrotherm plus (special voltage 100-127 V~)		502.0901.1
<p><b>ATMOS Aqua clean for rinsing and stimulation devices</b>  pressure-free water separation in accordance with EN 1717, water decontamination and discalcification  <i>Included in delivery:</i>  Device, connection material, 1 litre of ATMOS Green &amp; Clean WK disinfectant</p>		502.1200.0
Accessories		
<p><b>2 Water prefilter for ENT units, can be flushed out</b>  <i>Included in delivery:</i>  1 x reduction flange G3/4";  1 x nut - lip R1" T-chrome;  1 x solder fitting 90°;  1 x bracket for water filter;  2 x seal RD 17x24x2;  2 x seal RD 21x30x2.</p>		502.0990.0
<p><b>3</b> Rinsing connection for thermic stimulation, teflon application help for ear irrigation</p>		501.0355.0
<p><b>4</b> Jet connection, for water irrigation handle, L = 80 mm</p>		508.0427.0
<p>Jet connection, for water irrigation handle, L = 110 mm</p>		508.0429.0
<p><b>5</b> Splash protection, plexiglass</p>		501.0331.0
<p>Water supply hose G3/4", 1.5 m extension</p>		501.0315.6
<p>Waste water hose R 1/4", special length, surcharge/0.5 m</p>		502.0881.1
Consumables		
<p>Hose tips, 30 pcs. (single use),  minimum qty.: 2 packs (price per pack)</p>		502.0844.0
<p>O-ring for splash protection</p>		501.0331.1
<p>Filter cartridge for water filter</p>		502.0891.0
<p><b>6</b> Hygiene filter for water flushing  End positioned membrane filter for connection to the irrigation handle.  Medical product acc. to 93/42/EEC. Service life 4 weeks (Luer Lock).</p>		000.0918.0



1 ATMOS Hydrotherm plus




2 Water prefilter for ENT units



3 Rinsing connection for thermic stimulation



4 Jet connection for water irrigation handle



5 Splash protection



6 Hygiene filter for water flushing



**MedizinTechnik**

ATMOS MedizinTechnik GmbH & Co. KG  
Ludwig-Kegel-Str. 16  
79853 Lenzkirch / Germany  
Tel: +49 7653 689-370  
atmos@atmosmed.de

[www.atmosmed.com](http://www.atmosmed.com)