

ATMOS Variotherm plus

When precision counts!



ATMOS Variotherm plus

ATMOS Variotherm plus ear irrigation		REF
<p>1 Compact device for thermal nystagmus stimulation</p> <p>Thermal stimulation unit with electronically controlled flow heater, irrigation mode (37°C), cold- and warm-stimulation mode (each to be preset from 20° - 47° C) with time switching (1-99 s), standby function to reduce water and energy consumption , highest temperature accuracy and safety for the patient. Automatic flow limitation (150ml/min)</p> <p><i>Included in delivery:</i></p> <p><i>3 x nozzle attachment -II straight, 80 mm; 1 x splash protection; 1 seal for splash protection, 1 package of hose tips for nozzle, (30 pcs.); 1 x water filter housing R.V.-A1/4"; 1 x water filter CellPor-Pe; 1 x mains cable, 1x operating instructions.</i></p>		502.0900.0
ATMOS Variotherm plus (Sonderspannung 100-127 V~)		502.0900.1
<p>ATMOS Aqua clean for rinsing and stimulation devices</p> <p>Pressure-free water separation in accordance with EN 1717, water decontamination and discalcification</p> <p><i>Included in delivery:</i></p> <p><i>Device, connection material, 1 litre of ATMOS Green & Clean WK disinfectant</i></p>		502.1200.0
Accessories		
<p>2 Water prefilter for ENT units, can be flushed out</p> <p>Included in delivery:</p> <p>1 x reduction flange G3/4";</p> <p>1 x nut – lip R1" T-chrome;</p> <p>1 x solder fitting 90°;</p> <p>1 x bracket for water filter;</p> <p>2 x seal RD 17x24x2;</p> <p>2 x seal RD 21x30x2.</p>		502.0990.0
<p>3 Rinsing connection for thermic stimulation, teflon application help for ear irrigation</p>		501.0355.0
<p>4 Jet connection, for water irrigation handle, L = 80 mm</p>		508.0427.0
<p>Jet connection, for water irrigation handle, L = 110 mm</p>		508.0429.0
<p>5 Splash protection, plexiglass</p>		501.0331.0
<p>Water supply hose G3/4", 1.5 m extension</p>		501.0315.6
<p>Waste water hose R 1/4", special length, surcharge/0.5 m</p>		502.0881.1
Consumables		
<p>Hose tips, 30 pcs., single use)</p> <p>minimum qty.: 2 packs (per pack)</p>		502.0844.0
<p>O-ring for splash protection</p>		501.0331.1
<p>Filter cartridge for water filter</p>		502.0891.0
<p>6 Hygiene filter for water irrigation</p> <p>End positioned membrane filter for connection to the irrigation handle. Medical product acc. to 93/42/EEC.</p> <p>Service life 4 weeks (Luer Lock).</p>		000.0918.0



1 ATMOS Variotherm plus




2 Water prefilter for ENT units



3 Rinsing connection for thermic stimulation



4 Jet connection, for ear irrigation handle



5 Splash protection



6 Hygiene filter for water flushing



MedizinTechnik

ATMOS MedizinTechnik GmbH & Co. KG

Ludwig-Kegel-Str. 16

79853 Lenzkirch / Germany

Tel: +49 7653 689-370

atmos@atmosmed.de

www.atmosmed.com