

The perfect Full HD camera for office and outpatient clinic applications

ATMOS® CAM 41 HD

Single chip Full HD camera



The high sensitive CCD chip guarantees perfect images in combination with an ATMOS medical monitor, the patented ATMOS LED light sources, the ATMOS® Stroboscope and ATMOS® HD endoscopes.

Preinstalled profiles for different scopes and applications make working with this camera very simple and intuitive. Functions like image freeze, 180 degree image rotating or digital zoom help the doctor to concentrate on the essentials.

The ATMOS® CAM 41 HD fits perfectly into the well known ATMOS workflow and digital system integration, e.g. ATMOSoft 31/41, which captures pictures, movies and additional patient data.

Technical data	Resolution : Full HD 1080p / 1080i / 720p Automatic White Balance Electronic shutter : (1/50 to 1/100000) / 3 window sizes 2YC S- video output : (1 VCC & 0.3 VCC / 75) PAL or NTSC 2 Composite output : (1 VCC 75 oh) PAL OR NTSC DVI-D output, RGB output, Peripheral output - 2 foot switch input
-----------------------	---

Description	REF
ATMOS® CAM 41 HD incl. control unit, camera head and HD coupler	507.5500.0