


Professional face visors

# ATMOS Protection Shields



# Tailored Solutions for Medical and Nursing Staff



## The ATMOS Headpiece Professional — for maximum comfort and unrestricted mobility

The ATMOS Headpiece Professional has all the advantages of the standard version and offers added comfort and greater freedom of movement. Due to the special mounting system on this headpiece, the ATMOS Protection Shield Professional maintains an optimum distance to the face and can be flipped up if needed.

ATMOS Headpiece Professional + 100x ATMOS Protection Shields Professional  
REF: 507.4590.0

10x ATMOS Protection Shields Professional  
REF: 507.4591.0



## The ATMOS Headpiece — for medical and nursing staff

Together with the ATMOS Protection Shield, the individually adjustable headpiece offers the user simple but effective protection. The headpiece is reusable: the mounting strap is machine washable and the headpiece can be wipe disinfected.

ATMOS Headpiece + 100x ATMOS Protection Shields  
REF: 507.4585.0

10x ATMOS Protection Shields  
REF: 507.4586.0



### Notes:

The ATMOS Protection Shield is not an approved medical device. It is not guaranteed to block the transmission of viral pathogens, reduce the risk of disease infection or prevent any disease.

