


Tailored solutions
for medical applications

ATMOS

Protection Shields for devices



Three Tailored Solutions for Your Medical Applications



The Cordless ATMOS HL 21 LED Headlamp — with visor for expanding your ATMOS headlamp

With the ATMOS HL 21 LED, users have a cordless headlamp with supreme light intensity and a long battery life at their disposal. It can be expanded to include a transparent visor by means of the matching ATMOS Protection Shields. They can be easily secured to the headpiece of the headlamp.

ATMOS HL 21 LED Headlamp + 10x ATMOS Protection Shields
REF: 507.4505.0

10x ATMOS Protection Shields
REF: 507.4586.0



The ATMOS i View Protection Shield — for easy mounting on your ATMOS microscope

The transparent ATMOS i View Protection Shields can be mounted easily and quickly to the lens tube of the ATMOS i View. They thus form a barrier between the user and the patient without blocking the user's view.

10x ATMOS i View Protection Shields
REF: 538.4600.0

Note: Also fits OPMI pico



The ATMOS Endoscopy Protection Shield — attachment for your rigid endoscope

The transparent ATMOS Endoscopy Protection Shields can be attached to rigid endoscopes that have a diameter between 4 mm and 12 mm. During endoscopy, they thus form a barrier between the user and patient without blocking the user's view.

Also find out about the ATMOS iQam USB camera for displaying your endoscopic images digitally!

10x ATMOS Endoscopy Protection Shields
REF: 950.0270.0



Notes:

The ATMOS Protection Shield is not an approved medical device. It is not guaranteed to block the transmission of viral pathogens, reduce the risk of disease infection or prevent any disease.