



Easy to
use



Maximum
flexibility



Patented
folding

ATMOS Disposable Suction System

**Modern fluid
collection system
for any situation**

*Exactly
what
I need*

ATMOS Disposable Suction System



SAFE AND EASY HANDLING

- Easy removal of the full suction bag thanks to two handles
- Protection against oversuction and the ingress of bacteria thanks to the hydrophobic bacterial filter



MAXIMUM FLEXIBILITY

- Direct attachment to all electrical suction devices with an equipment mount holder according to DIN EN ISO 19054
- Attachment to any equipment rail using a rail clamp



PATENTED FOLDING OF THE BAG

- Up to 60% smaller packaging
- Lower transport costs for delivery and in-hospital transport
- Lower storage costs
- Significantly shorter priming time: easy insertion of the bag into the external canister



MedizinTechnik

Our ATMOS
product specialists
are available for
individual
consultation

ATMOS MedizinTechnik GmbH & Co. KG

Ludwig-Kegel-Str. 16
79853 Lenzkirch / Germany

Tel +49 (0)7653-689-0
Fax +49 (0)7653-689-190
E-Mail info@atmosmed.com

www.atmosmed.com

Your local contact: