



ATMOS iStore

The cabinet system for your examination room

*Exactly
what
I need*

iStore

intelligent | individual | infection control

HIGH INFECTION CONTROL STANDARDS

The choice of surfaces supports and facilitates compliance with infection control requirements and their implementation. The ATMOS iStore thus ensures that the guidelines are implemented excellently in practice operations.

Optimum
infection
control

PROMISING INTERIOR DESIGN

Our many years of experience in medical technology, the constant exchange of information with doctors, and a high degree of flexibility enable us to develop the best solutions for your practice.

In-depth
expertise

Individual
design

AN INTERPLAY OF DESIGN AND MATERIAL

Using only high-quality materials, the design is tailored to your needs. We make use of the available space in your practice rooms in the best possible way, integrating materials and technology that form a unity and save time.

Attractive
price

APPEALING PRICES FOR EVERY BUDGET

Thanks to the different materials and the high degree of modularity, individual solutions are also possible at a low price.




MULTITUDE OF DESIGN OPTIONS

Materials with ideal properties.

When selecting cabinet materials, quality is always our top priority. That's why we work with brand-name products that combine functionality, design, and durability. The materials we use meet the high requirements in the medical field and at the same time create a pleasant atmosphere for your patients and employees.

THE MATERIAL CORIAN®

CORIAN® is an innovative, durable, and high-quality mineral material. Thanks to its dirt- and bacteria-repellent properties, it is suitable for use in medical technology. The material has a nonporous surface and is especially hygienic. The work surfaces are extremely easy to clean and resistant to stains and chemicals. With our modules in the standard color Glacier White, we offer you a timeless design without compromise.



Surface of conventional materials, e.g. stainless steel CORIAN® Surface



DESIGN VERSATILITY

The rich palette of timeless colors and textures leaves nothing to be desired. Since CORIAN® can be shaped during the production phase, you can choose your own individual design. This makes it possible to create a uniform look, and the cabinet elements appear as one unit. With the decorative panels, we offer you an excellent alternative that likewise fulfills the required properties, even for low-cost concepts.



DECORATIVE PANEL MATERIAL

Our decorative panels are durable and hard-wearing. At the same time, they allow for a stylish and appealing design. Cleaning is easy, and you benefit from an inexpensive solution that meets all the requirements in your practice.



Individual design

info



To find detailed information on the mineral material CORIAN®, simply scan the QR code

HYGIENE-OPTIMIZED WORKFLOW

Unbeatable concept. Hygienically clean.


To minimize the risk of cross-contamination, we have developed an innovative concept for the correct handling of instruments. Clean and unclean instruments are strictly separated, and stringent infection control regulations and standards can be met at all times.



Optimum
infection
control

THE ATMOS INFECTION CONTROL CONCEPT FOR SPATIAL SEPARATION

You can take the disinfected and sterile instruments out of the cabinet and, once treatment is completed, dispose of them properly on the opposite side. The instruments are stored in airtight drawers equipped with Servo-Drive® electrical technology. This unique concept from ATMOS allows you to safely separate clean and unclean instruments.



DRAWERS FOR PERFECT INSTRUMENT MANAGEMENT

SERVO-DRIVE®

The drawers of the cabinets made of CORIAN® can be easily opened by simply tapping them with your knee, elbow, or back of your hand and closed again in the same way. This prevents the surfaces from being touched, thus virtually eliminating cross-contamination. Our cabinets made of decorative panels are available with the mechanical opening support called TIP-ON®, which allows you to open the drawers without handles by simply tapping them.



CLOSED SYSTEM

When equipped with the instrument and/or waste disposal, the cabinet can also be used to dispose of used instruments or other waste. The surface is made of CORIAN® and includes openings for easy disposal.



info Compatible with the following products from Aesculap (B. Braun)::

- Disposal container JK062R (150 mm stainless steel)
- Lid without opening JK068R
- Lid with opening JK067R

FURTHER OPTIONS



BACKSPASH / MILLED GROOVE


In addition to a flat worktop, our ATMOS iStore cabinet system offers you the option of choosing one with a backslash or milled groove. These features prevent liquids from running behind your cabinet system.

FINISHED SIDES

For the cabinet exteriors, the finished sides can be selected individually. This gives you limitless options when it comes to visual design. Determine the color and material of your cabinets' finished sides yourself and turn your visions into reality.

MOBILITY

If desired, you can configure the cabinet modules for mobile and room-independent use. With the help of practical swivel casters, our ATMOS iStore cabinets can be moved freely so that you can make the best possible use of your space.



Attractive
price

COST-EFFECTIVENESS AT ITS BEST

Excellent solutions. Fair prices.

Since quality, design and infection control are our focus, our product line creates the best possible solutions for any size room and budget. By combining the ATMOS iStore cabinets with our ENT treatment units, we offer you a modular product range that adapts to your wishes and financial possibilities.

PERFECT SOLUTIONS AT SMART PRICES

We provide you with state-of-the-art technology and an instrument management system that meets infection control requirements. We can offer all of this to you at an attractive price and therefore always find the best possible solutions for your rooms and needs. ATMOS iStore cabinets are thus a perfect complement to the ATMOS treatment units.



TWO MATERIALS FOR INFECTION CONTROL AND DESIGN

Our product lines focus on infection control, quality, and design. The material CORIAN® combines these three properties and is therefore a leader in the premium segment. With our line of cabinets made of decorative panels, we offer you a hygienic and visually appealing alternative that fits exactly into your budget.



INDIVIDUAL DESIGN OPTIONS

Thanks to the countless ways of combining our modules, you can arrange your work processes to make them even more efficient and economical. In addition, the various decors and design options offer you an exact custom fit while staying within your budget.



OUR STANDARD COLORS FOR DECORATIVE PANELS

We also offer you a large selection of decorative panels even if you decide on the lower-cost version.

FRONT DECORS



Alpine White
(W1100 ST9 NCS S0502/R50B)




Cream White
(W911 ST2 NCS S0502/G50Y)




Sand Beige
(U156 Texture 9 NCS 1505/Y50R)



Cashmere Grey
(U702 Texture 9 NCS S2005/Y60R)




Light Grey
(U708 ST9 NCS S200/N)



Grey
(U732 ST9 – Dust Grey, NCS S5000/N)



Tyrolean Blue
(U504 ST9)




Vicenza Oak
(H3157 ST12)



Cascina Pine
(H1401 ST22)



Deep Sea Blue
(U560 ST9)




Lincoln Walnut
(H1714 ST19)




Grey Corbridge Oak
(H3156 ST12)


WORKTOP DECORS



Pearl Grey
(U763 ST9)



Grey
(U732 ST9)



White Grey
(U775 ST9)

DECADES OF KNOW-HOW

Our expertise. Your solution.

For more than 100 years, we have been developing technical solutions for the medical field and supplying manufacturers, distributors, hospitals, and doctors worldwide with medical suction systems as well as complete solutions for ENT, gynecology, and swallowing diagnostics. The combination of medical and technical know-how enables us to respond exactly to your wishes and special requirements and to find suitable solutions for your everyday work.

OUR ROOTS IN BERLIN

ATMOS has been developing medical solutions since 1888. What began in a small pharmacy in Berlin has evolved over the years into an internationally active manufacturer of medical technology devices and products. The first patent for a self-developed medical device was granted as early as 1927. In the years that followed, medical research and development continued to advance. Thanks to our wealth of knowledge and experience, we are able to find the best technological solutions for your medical requirements.



EXPERTISE IN MEDICAL TECHNOLOGY AND CRAFTSMAN'S KNOW-HOW UNITED

In 2018, ATMOS MedizinTechnik expanded to include a cabinetry company. With ATMOS Interieur, we have now brought craftsmanship in-house. Innovative ideas, planning, and production are thus united under one roof, and your individual wishes can be implemented more quickly and effectively. As a result, we are able to fully equip all rooms in your practice.



WE TURN YOUR WISHES INTO SOLUTIONS

Thanks to many years of cooperation between our experts in the field of medical technology and the trades, our furniture program offers you solutions that perfectly complement and shape your practice and examination rooms. Our expertise and close cooperation with healthcare professionals are the basis for our work and the guarantee for your trust.



In-depth
expertise





ATMOS MedizinTechnik GmbH & Co. KG

Ludwig-Kegel-Str. 16
79853 Lenzkirch / Germany

Phone +49 (0)7653-689-0
Fax +49 (0)7653-689-190
E-mail info@atmosmed.com

www.atmosmed.com